



# MODEL 200

## Dimpling Grinder

An easy-to-use, state-of-the-art, mechanical thinning instrument designed for the reproducible preparation of high-quality electron microscopy specimens.

### Model 200 Dimpling Grinder Specifications

<b>Grinding control</b>	Automated grinding rate control Grinding force adjustable with micrometer counterweight system Independent control of grinding wheel and specimen rotation speeds
<b>Specimen stage</b>	Precise specimen stage rotation Magnetic mount allows easy specimen positioning
<b>User interface</b>	All program inputs via front panel keypad Specimen thickness indicated on a liquid crystal display
<b>Specimen illumination</b>	Specimen observation in either transmitted or reflected light
<b>Enclosure</b>	Weight: 18 lb (8.2 kg) Size: 8.2 in (208 mm) width x 6.5 in (165 mm) height x 13.5 in (343 mm) depth
<b>Power requirements</b>	110/220 V AC, 50/60 Hz, 375 W
<b>Warranty</b>	One year



E.A. Fischione Instruments, Inc.  
9003 Corporate Circle  
Export, PA 15632 USA  
Tel: +1 724.325.5444  
Fax: +1 724.325.5443  
info@fischione.com  
www.fischione.com

©2013 E.A. Fischione Instruments, Inc. All rights reserved.  
Document Number SP200 Revision 01 08/2013