



MODEL 160

Specimen Grinder

The Model 160 Specimen Grinder mechanically prethins specimens for transmission electron microscopy (TEM).

Model 160 Specimen Grinder Specifications

Specimen size	Up to 18 mm diam
Dimensions	3.0 in (76 mm) diam x 2.6 in (66 mm) height
Weight	2.4 lb (1.1 kg)
Warranty	One year



E.A. Fischione Instruments, Inc.
9003 Corporate Circle
Export, PA 15632 USA
Tel: +1 724.325.5444
Fax: +1 724.325.5443
info@fischione.com
www.fischione.com

©2013 E.A. Fischione Instruments, Inc. All rights reserved.
Document Number SP160 Revision 01 08/2013